

Blend-It™

Bench top sample blender

The Blend-It LC005 sample blender is ideal for blending composite samples prior to analysis.

The LC005 is not only rugged and simple to use but it will allow samples to be prepared under standard repeatable conditions.

The Blend-It can also be used for small scale manufacturing.



Blend-It LC005 with a 1 litre Mini-Drum

The adjustable clamp allows various container from 1 litre to 5 litres to be fitted

Blend-It LC005 features:

- Efficient asymmetric blending action
- Complete with adjustable drum clamp (for drums 1 to 5 litres)
- Speed fixed at 28 r.p.m.
- Timer - Range 0 to 30 minutes
- On/Off Switch
- Requires 230v 50Hz supply - (other voltages on request)
- Easily portable - weight 12kgs
- Carrying handles
- GMP correct design
- Comprehensive Operator and Maintenance Manual
- CE Mark
- All 304 grade stainless steel construction (contact parts 316 stainless steel)



1 litre and 2 litre Mini-drums for use with the LC005



Part No.	Description
8750L-01	Blend-It LC005 with drum clamp (clamp suitable for all mini-drums from 1 to 5 litres)
8750D-01	Mini-drum - 1 litre - 316L stainless steel - silicone 'O' ring seal - 3 toggle clamps
8750D-02	Mini-drum - 2 litres - 316L stainless steel - silicone 'O' ring seal - 3 toggle clamps
8750D-03	Mini-drum - 3 litres - 316L stainless steel - silicone 'O' ring seal - 3 toggle clamps
8750D-04	Mini-drum - 4 litres - 316L stainless steel - silicone 'O' ring seal - 3 toggle clamps
8750D-05	Mini-drum - 5 litres - 316L stainless steel - silicone 'O' ring seal - 3 toggle clamps
Please contact us for prices on other containers, options etc.	

Note: The volume being blended must not exceed two thirds of the total container volume

Larger Blenders

Sampling Systems' affiliate company Pharmatech Ltd can provide a range of other blenders right up to 6000 litre machines. Contact us for further information or visit:

www.pharmatech.co.uk

