

Chemical Dipper

Robust polypropylene Dipper for sampling aggressive liquids

Material: PP construction, ideal for aggressive chemicals

Notes:

1. Solid handle - non telescopic
2. Beaker is held securely - it can not be removed
3. Beaker has a pouring spout



Note: This sampler has crevices and so is not recommended for hygienic applications

Part No.	Nominal Beaker Volume	Overall Length
5065F-250	250ml	1000mm
5065F-500	500ml	1000mm
5065F-1000	1000ml	1000mm

NEXT DAY
DELIVERY

