

No Cleaning — No Cross Contamination — Reduce Costs

Powder Spatulas

Easy to use spatulas - made in a cleanroom

**ATEX Rated
Zones: 0,1 & 2**



Powder Spatulas are ideal for activities such as sampling and dispensing. They are manufactured from tough polystyrene and have a comfy ergonomic handle.



- Made from white Polystyrene (FDA grade - conforms to CFR 177.1640)
- Tough, efficient 'U' shaped blade
- Individually wrapped
- Manufactured and packed in a medical cleanroom (Class 100,000) to ensure cleanliness
- Two sizes
- ATEX rated for zones 0, 1 & 2
- Available pre-sterilised and non sterile
- Smooth ergonomic round handle



Powder Spatulas are robust and easy to use



Part No.	Blade Length	Overall Length	Pre-Sterilised	Quantity per Box
8095A-150	150mm	230mm	No	100
8095A-250	250mm	330mm	No	100
8095A-150S	150mm	230mm	Yes	100
8095A-250S	250mm	330mm	Yes	100

Discounts Available for Bulk Orders or Yearly Contracts