



**No Cleaning — No Cross Contamination — Reduce Costs**

# Thumb Tube Pipettes

## Low Cost Disposable Liquid Samplers

- Two styles available - Thumb Tube Pipettes & Jumbo Pipettes
- Simple to use
- Individually wrapped
- Assembled in a medical cleanroom
- Available pre-sterilised



### Thumb Tube Pipettes - small volume



Part No.	Overall Length	Overall Diameter	Internal Diameter	Pre-Sterilised	Material	Max Sample Volume*	Quantity per Box
5055A-50	1000mm	12mm	9mm	No	HDPE	60ml*	50
5055A-50S	1000mm	12mm	9mm	Yes	HDPE	60ml*	50

\* The maximum sample volume will depend upon how deeply the sampler is inserted into the product