

No Cleaning — **No Cross Contamination** — **Reduce Costs**

Disposable Dispensing Ware

Scoops with Lids



- Scoops made from rigid, white polystyrene
- Snap on lid to secure sample
- Pre-sterilised
- Individually wrapped



Part No.	Nominal Scoop Volume	Quantity per Box
8035A-60	60ml	100
8035A-120	120ml	100

ID-Scoops



- Made from white polyethylene
- Metallic strip in handle aids detection by metal detector
- Pre-sterilised
- Individually wrapped



Part No.	Nominal Scoop Volume	Quantity per Box
8030A-60	60ml	100

Sample Spoons



- Made from white polystyrene
- Individually wrapped in a re-sealable bag
- Pre-sterilised
- Cost effective



Part No.	Nominal Spoon Volume	Quantity per Box
8040A-120	3ml	100
8040A-150	10ml	100

Discounts Available for Bulk Orders or Yearly Contracts