

Easy Pumps

Economical and versatile pumps


The Easy Pump Range is ideal for many laboratory applications where small to medium quantities of liquid need to be pumped. Easy pumps are available either hand or foot operated.

Each pump is supplied with a Universal Bung - this allows the pump to be used with any container which has a round opening between 30 and 70mm.

These pumps are designed to be easy to install, easy to use and easy to clean.




Hand Operated Easy Pumps



Part No.	Tube Diameter	Maximum Length of Dip Tube	Tube Material	Capacity
5056A-18	18mm	700mm	Rigid PP	26L/min
5056A-15	15mm	700mm	Rigid PP	20L/min
5056A-12	12mm	700mm	Rigid PP	14L/min
5056B-16	16mm	700mm	Rigid PVDF	20L/min
5056C-18	18mm	900mm	304 Stainless steel	26L/min
5056D-10	10mm	700mm	Flexible PVC	12L/min
5056E-10	10mm	450mm	Rigid PP	4L/min

Every EasiPump is supplied with a Universal Bung - for use with openings 30 to 70mm in diameter

Foot Operated Easy Pumps



Part No.	Tube Diameter	Maximum Length of Dip Tube	Dip Tube & Delivery Tube Material	Capacity
5058A-22	22mm	900mm	Rigid PP	30L/min
5058A-18	18mm	700mm	Rigid PP	26L/min
5058A-15	15mm	700mm	Rigid PP	20L/min
5058B-16	16mm	700mm	Rigid PVDF	20L/min
5058C-18	18mm	900mm	304 Stainless steel	26L/min
5058D-10	10mm	700mm	Flexible PVC	12L/min

Every EasiPump is supplied with a Universal Bung - for use with openings 30 to 70mm in diameter