

# NIR


## PharmaView

### Real-time Moisture Sampler for Fluid Bed Driers

- Real-time moisture measurements
- Non-invasive sampling technique
- Can be retro-fitted to an existing Fluid Bed Drier
- Uses Near Infra Red measurement system
- Self cleaning lens for consistent performance
- Unique patented design
- GMP correct design and construction
- **Never over-dry a batch again**



Cross section through a FBD bowl showing the PharmaView



### Lens Cleaning Systems

The PharmaView has a unique Lens Cleaning System that allows the device to ensure all debris is removed from the lens before a reading is taken.

## Call now for further information