

PTFE Dipper

Highly resistant Liquid Sampler

The PTFE Dipper is Ideal for sampling aggressive chemicals.

Suitability:	Liquids (low, medium & high viscosity)
Material:	Sampling Cup: PTFE Handles & extension rods: PTFE coated steel rod
Operation:	Insert Dipper into liquid. Take sample



To increase the effective length of the sampler use extension rods.

Note: Extension rods can only be used with the 600mm long removable handles.



Part No.	Nominal Volume	Cup Diameter	Cup Height	Handle
5100A-10	10ml	25mm	34mm	with 155mm long handle - fixed in place
5100A-100	100ml	55mm	65mm	600mm long - removable
5100A-250	250ml	70mm	95mm	600mm long - removable
5100A-500	500ml	80mm	125mm	600mm long - removable
5100A-1000	1000ml	105mm	155mm	600mm long - removable
5100E-600	Extension rod - 600mm long (not suitable for use with 5100A-10)			