

# Sieve-It™

Robust Sieve Shakers for particle size analysis



## Sieve-It™ Octagon

The Sieve-It Octagon is a high performance test sieve shaker offering excellent operator control for maximum efficiency.

- Efficient, reproducible sieving action
- Robust compact design
- Quick release clamping system
- Power supply: 230v 50Hz or 110v 60Hz  
Other voltages & frequencies on request
- Digital timer. Vibration amplitude control
- Continuous or intermittent operation
- Minimal maintenance
- Low noise during operation. Anti-vibration mountings
- Weight: 43kg
- Accommodates up to 8 full height sieves
- CE marked



Part No.	Description
8700A-01	Sieve-It Octagon - 230v 50hz
8700A-02	Sieve-It Octagon - 110v 60Hz



## Sieve-It™ Minor

The Sieve It Minor is ideal for low - medium use.

- Efficient, reproducible sieving action
- Compact design
- Quick release clamping system
- Power supply: 230v 50Hz or 110v 60Hz  
Other voltages & frequencies on request
- Simple timer. Vibration amplitude control
- Continuous or intermittent operation
- Minimal maintenance
- Low noise during operation. Affordable
- Light Weight only 17kg
- Accommodates 8 off sieves
- CE marked



Part No.	Description
8700B-01	Sieve-It Minor - 240v 50hz
8700B-02	Sieve-It Minor - 110v 60Hz



## Test Sieves

High Quality test sieves for accurate and efficient particle size analysis.

- Woven wire mesh sieves
- Construction: All stainless steel or brass frame and stainless steel mesh
- Range of aperture sizes available (4mm to 20 microns)
- Certificate of Conformity with each sieve
- Other styles of sieve available on request (Perforated Plate Sieves, Microplate Sieves, Grain sieves etc)
- Test sieves manufactured to every major International standard (inc. ISO, BSI, & ASTM)

See next page for sizes

# Sieve It™ cont.

## Test Certificates

Every sieve is supplied with a Certificate of Conformity.



## Sieve Height

The standard height for a 200mm diameter sieve is 50mm. Half height sieves are available - ask for details.

Other diameter sieves have different heights - ask for details



### 200mm Diameter Test Sieves

All Stainless Steel Sieves		Brass Frame & Stainless Mesh Sieves	
Part No.	Aperture	Part No.	Aperture
8700S-01	4mm	8700B-01	4mm
8700S-02	2mm	8700B-02	2mm
8700S-03	1mm	8700B-03	1mm
8700S-04	850 microns	8700B-04	850 microns
8700S-05	700 microns	8700B-05	700 microns
8700S-06	600 microns	8700B-06	600 microns
8700S-07	500 microns	8700B-07	500 microns
8700S-08	425 microns	8700B-08	425 microns
8700S-09	300 microns	8700B-09	300 microns
8700S-10	250 microns	8700B-10	250 microns
8700S-11	180 microns	8700B-11	180 microns
8700S-12	150 microns	8700B-12	150 microns
8700S-13	125 microns	8700B-13	125 microns
8700S-14	106 microns	8700B-14	106 microns
8700S-15	90 microns	8700B-15	90 microns
8700S-16	75 microns	8700B-16	75 microns
8700S-17	53 microns	8700B-17	53 microns
8700S-18	32 microns	8700B-18	32 microns
8700S-19	25 microns	8700B-19	25 microns
8700S-20	20 microns	8700B-20	20 microns

Other sieve diameters available on request