

Stainless Dipper

GMP correct Liquid Sampler

A high quality, crevice-free dipper suitable for sampling a wide range of liquids.

Suitability: Liquids (low, medium & high viscosity)

Material: 316L stainless steel

Operation: Insert Dipper into liquid. Take sample



Part No.	Nominal Volume	Overall Length	Outside Cup Diameter
5005B-100	100ml	1000mm	50mm
5005B-250	250ml	1000mm	65mm
5005B-500	500ml	1000mm	90mm

Other handle lengths available on request