



Mini-Pots

- Multi purpose small volume containers
- 316 stainless steel
- Crevice free interior
- Screw on lid
- Satin finish
- Water tight - no extra seal required

Part No.	Volume	Overall Diameter	Height
A426-25	25ml	38mm	45mm
A426-50	50ml	50mm	45mm
A426-100	100ml	50mm	71mm
A426-250	250ml	75mm	75mm
A426-500	500ml	100mm	75mm
A426-750	750ml	100mm	115mm
A426-1000	1000ml	100mm	155mm



Pellet Pots

- Multi purpose small volume containers
- 316 stainless steel
- Crevice free interior
- Screw on lid
- Satin finish

Part No.	Volume	Overall Diameter	Height
A056-7	7ml	16mm	41mm

Special Sizes on Request -
ask for details



Earthing lug welded to the handle of a 500ml PharmaScoop

Earthing Lugs

- Extra cost option for any stainless steel sampler or item of dispensing ware
- Provides a convenient connection point for earthing (grounding) wires
- Ideal when working with flammable solvents or in ATEX rated environments
- Earthing lugs can only be welded to almost any stainless steel device

Part No.	Description
LUG-1	Earthing Lug - 316L stainless steel. Includes welding to device

Mini-Drums & Storage Containers

- Multi purpose containers
- Crevice free interior
- Close fitting recessed lid with handle
- Flat bottom with chime ring to give added protection
- Complete with lid (toggle clamps & 'O' ring are extra)

Handles

- MiniDrum - Swing handle
- Storage Container - Two sturdy handles fitted on patch plates

Options

Toggle Clamps -



To secure the lid we recommend the use of toggle clamps

Use 3 toggle clamps per Mini-Drum

'O' Ring Seals -



Available in:
Silicone
Viton
Nitrile



304 Stainless Steel Part No.	316L Stainless Steel Part No.	Nominal Volume	Height	Diameter
A434-1	A446-1	1 litre	155mm	100mm
A434-2	A446-2	2 litres	175mm	130mm
A434-3	A446-3	3 litres	185mm	155mm
A434-5	A446-5	5 litres	240mm	175mm
A434-6	A446-6	6 litres	280mm	175mm
A434-8	A446-8	8 litres	285mm	200mm
A434-10	A446-10	10 litres	295mm	220mm
A434-12	A446-12	12.5 litres	285mm	240mm
A434-15	A446-15	15 litres	285mm	270mm



304 Stainless Steel Part No.	316 Stainless Steel Part No.	Nominal Volume	Height	Diameter
A454-20	A466-20	20 litres	387mm	270mm
A454-25	A466-25	25 litres	387mm	300mm
A454-30	A466-30	30 litres	455mm	300mm
A454-50	A466-50	50 litres	555mm	350mm
A454-75	A466-75	75 litres	620mm	400mm
A454-100	A466-100	100 litres	670mm	450mm



Pallet Knives

- All stainless steel construction
- Crevice free
- Satin finish

Part No.	Blade Length	Blade Width
A330-100	100mm	20mm
A330-150	150mm	25mm
A330-200	200mm	32mm
A330-250	250mm	38mm
A330-300	300mm	45mm
A330-350	350mm	



Scrapers

- All stainless steel construction
- Crevice free
- Satin finish

Part No.	Blade Width	Overall Length
A350-75	75mm	232mm
A350-100	100mm	275mm



Sieving Tools

- Ideal for pushing powders through mesh screens
- 316L stainless steel
- Radiused corners
- Crevice free construction
- Satin finish
- Robust handle

Part No.	Disc Diameter	Handle Length
A516-120	120mm	240mm



Stainless Spoons

- Ideal for transferring small amounts of powder
- Crevice free design
- Bright polished finish
- 304 stainless steel

Part No.	Nominal Volume	Length
A504-1	2ml	100mm
A504-2	10ml	150mm
A504-3	2ml	210mm
A504-4	10ml	300mm
A504-5	10ml	500mm



Powder Spoons

- Ideal for transferring small amounts of powder
- Crevice free design
- Bright polished finish
- 304 stainless steel

Part No.	Nominal Volume	Length
A100-20	20ml	170mm



Buckets

- 316L or 304 stainless steel
- Crevice free design
- Graduated
- Swing handle

Part No.	Nominal Volume	Stainless Steel
A264-6	6 litres	304
A264-8	8 litres	304
A264-10	10 litres	304
A264-12	12.5 litres	304
A264-15	15 litres	304
A276-10	10 litres	316L
A276-12	12.5 litres	316L
A276-15	15 litres	316L



Shallow Trays

- 304 or 316L stainless steel
- Crevice free construction
- 10mm deep

Part No.	Stainless Steel	Length	Width
A224-220	304	220mm	150mm
A224-290	304	290mm	210mm
A224-435	304	435mm	275mm
A224-520	304	520mm	275mm
A224-520A	304	520mm	350mm
A226-220	316L	220mm	150mm
A226-290	316L	290mm	210mm
A226-415	316L	415mm	240mm
A226-520	316L	520mm	350mm



Deep Trays

- 304 or 316L stainless steel
- Crevice free construction

Part No.	Stainless Steel	Dimensions (L x W x D)
A234-220	304	220 x 150 x 45mm
A234-355	304	355 x 250 x 55mm
A234-395	304	395 x 280 x 70mm
A234-475	304	475 x 315 x 90mm
A234-530	304	530 x 325 x 100mm
A236-220	316L	220 x 150 x 45mm
A236-355	316L	355 x 250 x 55mm
A236-395	316L	395 x 280 x 70mm
A236-475	316L	475 x 315 x 90mm



Buckets with Spouts

- 304 stainless steel
- Crevice free & satin finish
- Robust construction

Part No.	Volume	Stainless Steel
A404-10	10 litres	304
A404-15	15 litres	304
A404-20	20 litres	304

Note: Bucket shown above has an optional lid - ask for details



Square Scoops

- Heavy duty construction
- 316 stainless steel
- Crevice free
- Radiused corners for easy cleaning

Part No.	Nominal Volume
A366-1	2.5 litres

PharmaScoops



- Heavy duty construction
- 316 stainless steel
- Crevice free
- Satin finish
- Radiused corners for easy cleaning
- Huge range of sizes

Part No.	Nominal Capacity	Body Diameter	Body Length
A316-10	10ml	25mm	45mm
A316-50	50ml	38mm	70mm
A316-100	100ml	50mm	100mm
A316-500	500ml	95mm	140mm
A316-750	750ml	95mm	210mm
A316-1000	1000ml	110mm	235mm
A316-1250	1250ml	125mm	195mm
A316-1700	1700ml	135mm	230mm
A316-2500	2500ml	185mm	215mm



Ladles

- All stainless steel construction
- Crevice free
- Satin finish
- 316 stainless steel for 50 and 150ml size

Part No.	Volume	Material
A324-50	50ml	304 ss
A324-120	120ml	304 ss
A324-150	150ml	304 ss
A324-200	200ml	304 ss
A326-50	50ml	316 ss
A326-150	150ml	316 ss



Bailers

- 316 stainless steel construction
- Crevice free
- Satin finish

Part No.	Nominal Volume	Bowl Diameter
A296-1	1 litre	150mm
A296-2	2 litres	200mm
A296-4	4 litres	250mm



Economy Scoops

- 304 stainless steel
- Light duty construction

Note: The Economy Scoop is NOT suitable for GMP applications as there are crevices

Part No.	Nominal Capacity	Diameter
A304-120	0.12 litres	50mm
A304-450	0.45 litres	75mm
A304-1000	1 litre	105mm
A304-1500	1.5 litres	120mm
A304-2000	2 litres	150mm



Funnels

- 304 or 316 stainless steel
- Crevice free construction
- Fluted spout to aid product flow

Part No.	Stainless Steel	Top Diameter	Spout Diameter
A244-50	304	50mm	6mm
A244-100	304	100mm	16mm
A244-120	304	120mm	16mm
A244-150	304	150mm	18mm
A244-200	304	200mm	18mm
A244-250	304	250mm	18mm
A244-310	304	310mm	18mm
A246-150	316	150mm	18mm



Narrow Top Jugs

- 304 or 316 stainless steel
- Crevice free construction
- Non graduated except for 5L size

Part No.	Stainless Steel	Nominal Volume
A214-600	304	600ml
A214-1000	304	1000ml
A214-3000	304	3000ml
A214-5000	304	5000ml
A214-5000G	304	5000ml
A216-3000	316L	3000ml
A216-5000G	316L	5000ml

G = Graduated



Paddles

- 316 stainless steel
- Crevice free & satin finish
- Robust construction

Part No.	Overall Length	Blade Width	Blade Length
A344-1000	1000mm	125mm	240mm
A344-1500	1500mm	150mm	300mm